


최고의 know-how와 기술! Needle 전문기업!!



• Special design for effectiveness of sterilization.

• Design of ergonomical cap.

• Electric discharge of all products on hub hole.




Bevel Indicating mark

• The highest quality of silicone coating substantially reduces pain and smooth repeated injection


• Tri-bevel point, perfected treatment of surface protect the gum without damage of tissue

- 완벽한 표면처리, 잇몸조직의 손상과 통증을 최소화하는 **이상적인 연막각도!**
- 국내 및 해외 최초! **수출마취 전용 30G X-Long(30mm)** 개발!
- 바늘표면의 **실리코닝과 특수열처리**로 수 차례의 Injection에도 부드럽게 삽입!
- 총생산량의 95%를 수출하고 있으며 미국, 영국, 이태리 등 유럽을 비롯한 **전세계 100여 개국에서 최고의 품질로 인정**
- **국내 최초 31G(0.25mm)**제품 식약처허가 및 양산!

C-K MEDICAL / DENTAL NEEDLE POINT




BETWEEN BEVEL SURFACES



- Design of ergonomical cap.
- Special design for effectiveness of sterilization.
- Electric discharge of all products on hub hole.
- Tri-bevel point for easy tissue penetration without coring.
- The round cutting edge protecting the gum
- The most strict grinding angle, ratio of tri-bevel points reduce pain
- The highest quality of silicone coating substantially reduces pain and smooth repeated injection.
- Bevel Indicating mark

DENTAL INJECTION NEEDLES



SIZE	REF No. (INCH : I. METRIC : M)
27G (0.4)	21mm 27-21-(M)
	25mm 27-25-(M)
	30mm 27-30-(M)
	35mm 27-35-(M)
	38mm 27-38-(M)
30G (0.3)	40mm 27-40-(M)
	12mm 30-12-(M)
	16mm 30-16-(M)
	21mm 30-21-(M)
	25mm 30-25-(M)
31G (0.25)	30mm 30-30-(M)
	12mm 31-12-(M)

